Cowslips



Leading prosthetic device for cows that allows relief of the damaged claw

Features:

- The glue is mixed inside the shoe for simple, clean and easy application.
- Increases the mobility of the affected animal and allows for increase in milk production and weight gain.
- Can be easily cut and adapted for different shapes of feet and treatments.
- The unique form assures good surface contact with the claw allowing it to stay on for 6-8 weeks before shedding. naturally.





Cowslips

Product Information:

Cowslips are one of the leading prosthetic devices for lameness in Cows. Cowslips were developed in the UK by Dr Ernie Logan, and are now used worldwide.

Cowslips are fitted to the unaffected claw. First, the affected hoof should be lifted and cleaned, and the cause of lameness determined. If the unaffected claw is healthy, the cowslip can be applied to this claw – this allows load bearing relief on the damaged claw so healing can occur. The Cowslip will slough off of its own

accord in six to eight weeks, by which time the affected claw should have recovered. Use the Cowslip 'Original' on regular sized cows, 'Plus' for the large and heavy cows and 'X-large' for very large and heavy animals. Refer to the table below for further size information.

Code	Product
211583	Cowslips Original Mixed L+R 10-pack
211584	Original Left only 10-pack
211585	Original Right only 10-pack
211589	Cowslips Plus Mixed L+R 10-pack
211590	Cowslips Plus Left only 10-pack
211591	Cowslips Plus Right only 10-pack
211592	Cowslips X-Large Mixed L+R 4-pack

Code	Product
211600	Cowslips Powder Single App'n each
211601	Cowslips Liquid Single App'n each

ORIGINAL 116mm long For use on Jerseys, • 60mm wide Cowsli Guernseys, Heifers and • 16mm high Ayrshires. 129mm long (13mm longer than the original) PLUS 61mm wide For use on Holsteins, (1mm wider than the original) Friesians and Beef Breeds. 17mm high (1mm higher than the original) Cowslips

XL For use on large/abnormal claws, North American Holsteins. • 142mm long (26mm longer than the original) • 62mm wide (2mm wider than the original)

Controlic Ltd.

201600

• 18mm high (2mm higher than the original)



201601

www.shoofvet.co.nz New Zealand 0800 800 801