

Speedy Flexi Probe



Speedy Feeders

Features:

- Control the animal and the fluid delivery
- Bottle can be placed on ground
- Convert your speedy bottle to a drencher

Product Information:

The Speedy Flexi Probe allows for easy positioning of the probe without encumbrance of the bottle which can now be placed on the ground.

Developed to be compatible with Shoofs famous 2.5L & 4L Speedy bottles. Provides the user with maximum control over the animal and over

the fluid delivery by ensuring the appropriate fluid volume is administered, and at the correct adjustable flow rate (3-speeds) and cut-off as required. Simply change the speedy lid assembly from teat to flexi probe and convert your bottle feeder into a drencher.

The smooth clear flexible tubing is 110cm long with the last 28cm rigid with a specially designed probe tip reduces the possibility of malpositioning to the lungs. Colostrum is notorious for being thick, sticky and difficult to clean - the wide mouth of the speedy feeder bottle ensures quick and easy cleaning. An extremely versatile addition to your calf rearing kit.

Code	Product
224 582	Calf Drencher Speedy Flexi Probe Assy
207 404	Calf Feeder Speedy Feeder 2.5L cpt
213 731	Calf Feeder Speedy Feeder 4L cpt